


3L Cell Culture Bioreactor Systems, Unjacketed

Agitator Assembly for
Use with

CLS-1399-100
Overhead Stirrer.

Additional Adapters
are available for other
style motors.



18 Port Headplate (4)

Graduated Vessel (1)

Support Stand (2)

Vortexing Impeller (3)

Cell culture bioreactor systems include all the necessary base components needed in one simple package, supplied with or without a motor. By using motor adapters, the system can be retrofitted to work with almost any motor. These systems can be purchased as configured or customized to fit your specific application.